

Adenosine

Indication: SVT
Concentration = 3 mg/ml

COLOR	Doses (mg)	Give (ml)
Gray	0.45 mg	1 st 0.15 ml
	0.9 mg	2 nd 0.3 ml
Pink	0.66 mg	1 st 0.22 ml
	1.35 mg	2 nd 0.45 ml
Red	0.9 mg	1 st 0.3 ml
	1.8 mg	2 nd 0.6 ml
Purple	1 mg	1 st 0.33 ml
	2mg	2 nd 0.67 ml
Yellow	1.35 mg	1 st 0.45 ml
	2.7 mg	2 nd 0.9 ml
White	1.7 mg	1 st 0.6 ml
	3.4 mg	2 nd 1.2 ml
Blue	2.1 mg	1 st 0.7 ml
	4.2 mg	2 nd 1.4 ml
Orange	2.7 mg	1 st 0.9 ml
	5.4 mg	2 nd 1.8 ml
Green	3.3 mg	1 st 1.1 ml
	6.6 mg	2 nd 2.2 ml
40kg	4mg	1 st 1.3 ml
	8mg	2 nd 2.7 ml
45kg	4.5 mg	1 st 1.5 ml
	9 mg	2 nd 3 ml

Note: Immediately follow with a rapid 10-20ml NS bolus.

